

DOC Cases with 2 accounts

With DOC cases you may sometimes notice that you see priors listed on the patient history page, but even when you have all priors selected, the thumbnails do not show. Usually the reason is that the prisoner has been transferred so when the images are processed it adds the initials of the new prison behind the MRN so it doesn't match the original MRN. In these instances you will notice that the patient is listed again in the priors list rather than one long list of priors like the screen shot below.

The screenshot shows a medical interface with the following details:

- Top Left:** Patient information including "DOB: November 22, 1983", "PID: 1234CC", "Note: None", and "Alert: None". A red arrow points to a label "Same Name".
- Top Right:** Study title "X-RAY CHEST 2 VIEWS", "Routine" dropdown, "June 03, 2014", "Acc #:" (blacked out), and "Status: Report Available - ZZ".
- Below Top:** "Auto-Relevancy: Same Body Part" and "Set as Default".
- Table:** A table with columns: Rel, Study Date, Modality, Study Description, State, and Normalized Code. The first row shows "09:10", "DX", "X-RAY CHEST 2 VIEWS", "71020", and "95%".
- Below Table:** Patient information including "DOB: November 22, 1983 (30Y)", "ORC: SCC", and "1234BC". A red arrow points to "1234BC". A "Same Patient" checkbox is present.
- Table Rows:** Three rows of study data:
 - May 06 - 08:46 DX X-RAY CHEST 2 VIEWS 71020
 - Apr 08 - 11:30 DX X-RAY CHEST 2 VIEWS 71020
 - Mar 13 - 09:01 DX X-RAY CHEST MINIMUM 4... 71030

You will also see a percentage. This represents what the likelihood is that it is the same patient. If it is 95% it is most likely the same patient.

This close-up shows the "Normalized Code" column with the value "71020" and a "95%" likelihood percentage. Below the percentage is a "Same Patient" checkbox.

Look at the name and DOB, also for DOC patients the MRN will be the same except for the prison initials following the numbers.

Same.Name ←
 DOB: November 22, 1983
 PID: **1234CC**
 Note: None
 Alert: None

X-RAY CHEST 2 VIEWS
 Routine June 03, 2014
 Acc #: ██████████ Status: Report Available - ZZ

Auto-Relevancy: Same Body Part Set as Default

Reset Filters

Rel	Study Date	Modality	Study Description	State	Normalized Code
Filter...	Filter...	Filter...	Filter...		
	09:10	DX	X-RAY CHEST 2 VIEWS	<input checked="" type="checkbox"/>	71020
					95%
				<input type="checkbox"/>	<input type="checkbox"/> Same Patient
<input type="checkbox"/>	May 06 - 08:46	DX	X-RAY CHEST 2 VIEWS	<input checked="" type="checkbox"/>	71020
<input type="checkbox"/>	Apr 08 - 11:30	DX	X-RAY CHEST 2 VIEWS	<input checked="" type="checkbox"/>	71020
<input type="checkbox"/>	Mar 13 - 09:01	DX	X-RAY CHEST MINIMUM 4...	<input checked="" type="checkbox"/>	71030

Same, Name ←
 DOB: November 22, 1983 (30Y)
 PID: **1234BC**
 ORG: SCC

If you click on the box next to same patient the thumbnails of those priors will then become available to view as priors.

Same.Name PID: **1234BC** DOB: November 22, 1983 (30Y) ORG: SCC

Same Patient

<input checked="" type="checkbox"/>	May 06 - 08:46	DX	X-RAY CHEST 2 VIEWS	<input checked="" type="checkbox"/>	71020
<input checked="" type="checkbox"/>	Apr 08 - 11:30	DX	X-RAY CHEST 2 VIEWS	<input checked="" type="checkbox"/>	71020
<input checked="" type="checkbox"/>	Mar 13 - 09:01	DX	X-RAY CHEST MINIMUM 4...	<input checked="" type="checkbox"/>	71030

Relevant: Jun 3, 2014 X-ray-X-RAY CHES... May 6, 2014 X-ray-X-RAY CH... Apr 8, 2014 X-ray-X-RAY CH... Mar 13, 2014 X-ray-X-RAY CH...

Jun 3, 2014 X-ray - X-RAY CHEST 2 VIEWS

Same.Name
 DOB: Nov 22, 1983
 PID: **1234cc**
 ACC: ██████████
 VIEWED

We only merge these patient if the DOC asks us to.